



Adaptors and Drain Valves

Urocare® Quick Drain Valve™

Only a limited amount of finger or hand control is required to operate these valves independently. Finger loops or a crossbar may be attached to the lever to suit individual requirements—see *Figure A*. Fits Urocare® Reusable Latex Leg Bags and Urocare® Male Urinal Sheaths only.

Made with pliable, white-rubber tubing and plastic—**no metal parts to corrode**. May be used as a replacement for the Standard Bottom Drain Valve (#6004) preassembled on all Urocare® Reusable Latex Leg Bags and Urocare® Male Urinal Sheaths.

KEY FEATURES:

- Lever operated for quick, controlled drainage
- Threads are sealed to prevent urine from touching the skin while draining the bag
- All working parts are urine free
- No metal parts or sharp corners
- Tube may be cut to desired length



FIGURE A - "T"-Bar

An effective crossbar can be made by cutting 2" (5cm) from the bottom of the drain valve tube and pushing one end through the hole in the lever until it is centered, forming a "T" shaped bar.

DEHP-FREE

LATEX-FREE



LEARN MORE



Sold by the each, 1 valve per package.
These are typically stocked items.
Non-Sterile, DEHP & Latex-Free.

Item No.	Description	Size†	Case Item No.
700172	Extra-Large	0.5" x 72" (183cm) Tube Length	Not Available
700118	Large	0.5" x 18" (45cm) Tube Length	Not Available
700110	Standard	0.5" x 10" (25cm) Tube Length	7011012 (12/CS) 700110200 (200/CS)

† Approximate dimensions in inches(in) and centimeters(cm).