

UROCARE



ADHESIVE REMOVER PAD

CITRUS SCENTED WIPE

Safe • No Harsh Solvents • Non-Toxic
Biodegradable • Non-Irritating
No Artificial Dyes, Colors or Fragrance

Removes most types of adhesive from
skin and appliances
with Aloe Vera

For External Use Only
Keep out of reach of children

ADHESIVE REMOVER PAD #5600

DIRECTIONS: Apply directly to edge of adhesive while gently pulling on adhesive strip or appliance until separated from skin. Wipe off any remaining residue. Rinse with soap and water.

USES:

- Adhesive Tapes
- Adhesive Residue
- Ostomy Appliances
- Tar Removal
- TENS Accessories
- Nitro Patches

Questions/Comments? 1-800-423-4441

CONTAINS: Natural Oils, Aloe Vera, Lanolin and Organic Wetting Agents. Natural Citrus Scent.

UROCARE PRODUCTS, INC.
Pomona, CA 91767, U.S.A.



www.urocare.com